

Abstract

Time-Domain nuclear magnetic resonance (TD-NMR) represents an attractive alternative method for analyzing petroleum products such as motor oil and diesel fuel. This is due to its ability to analyze samples with little or no sample preparation, allowing fast data collection, and being a nondestructive and noninvasive method. Current TD-NMR methods that have been developed to analyze petroleum product samples utilize very large and expensive TD-NMR benchtop equipment. WaveGuide is developing a portable TD-μNMR (Time-Domain micro-NMR) that can be utilized in the field to determine the viscosity and authenticity of branded petroleum products.

TD-NMR was explored as a rapid method for simultaneous assessment of the quality parameters in conventional and synthetic motor oil samples. Data obtained with the relaxation decay curves employing a Carr-Purcell-Meiboom-Gill (CPMG) pulse sequence revealed tight and well-separated clusters of the motor oils, allowing discrimination of the motor oil samples according to their viscosity content and brand: 0W-30, 5W-30, 10W-30 and 10W-40.¹ A set of 10 ASI standards for sulfur in diesel fuel was also analyzed, using a CPMG pulse sequence.² The resulting data also showed well-separated clusters, allowing discrimination among the standards with sulfur weight percentage (wt %): 0, 0.1, 0.25, 0.5, 1, 1.5, 2, 3, 4, 5.

WaveGuide's portable TD-μNMR is the first battery powered instrument that is fully automated and robust enough to operate in the field without the need for a trained NMR technician. The device is 22x smaller, ~80x lighter, and 30x lower cost than contemporary commercial analytical systems.



How Does This Fit with PANIC's Mission?

WaveGuide's portable TD-NMR (time domain micro-NMR) provides solutions for field testing of petroleum products for counterfeiting. Testing can now be conducted on site at ports of entry where counterfeit products are likely to enter the supply chain. WaveGuide's TD-NMR allows non-scientists to conduct field testing of counterfeit motor oils as well as high sulfur contaminated diesel products. WaveGuide is providing a practical industry-focused solution by bringing TD-NMR to the masses.

Introduction

On a global scale, the counterfeiting business is highly profitable across multiple retail markets. Retail goods that are routinely counterfeited include high end clothes, jewelry, and cell phones. Counterfeiters create these fake replicas to answer the growing consumer demand for these products. The business of creating fake goods is not just for the luxury items as it extends into the energy industry, specifically with motor oil and other automotive lubricants. Operations of massive oil counterfeiting can be seen within the United States as well as countries like Russia, where in December of 2015, a large counterfeit oil ring was discovered outside of Moscow by Russian law enforcement, collecting a profit of \$164 million a year.^{3,4} Large brand names became damaged with the illegal activity such as Mobil, Castrol, Shell, and others.⁵ The individuals behind the ring even faked OEM oil with major car brands such as Toyota, Ford, and Mazda, selling the inauthentic oil under their brand names. The frequency of such counterfeits is both detrimental to the brands' reputations but also could pose potentially catastrophic and deadly consequences to the consumers.⁶ The risk of damage can be minimized if products can be analyzed and authenticated prior to falling into the customers hands. A device that can test a sample on the spot and within minutes determine if the product is genuine or counterfeit can save customers, vendors, corporations, and law enforcement significant time and money typically used to derail such counterfeit operations.⁴ The WaveGuide Formula™ is a truly portable, handheld device that can be used by nearly anyone.



The WaveGuide Formula™ portable μNMR
Little to no sample preparation required
Sample analysis take only minutes
Less than 50μL of sample typically used per test

Materials & Methods

Objective 1: Analyze a set of 19 commercially available motor oil samples (Table 1) and explore correlation between T2 and reported viscosity to correctly identify commercial samples.

Experiment Methods: Each sample was freshly prepared and tested (using 40μL sample volume) using the protocol: 3 runs each sample, accumulating 4 scans per run, scan duration of 0.4 seconds, a recycle delay of 2 seconds, and 400 μs for the echo period in the CPMG⁴ pulse sequence at 25°C. The magnetic field of the μNMR was at 0.50 Tesla. The data was fitted using a bi-exponential algorithm.

Objective 2 Analyze a set of crude oil samples (Table 3) to correlate T2 with the determined wt% asphaltene which is a crude oil contaminant that usually requires a 24-hour labor intensive method to determine to create a fast replacement test.

Objective 3: Analyze a set of crude oil samples (Table 2) with determined viscosity and density to provide a faster analytical method for determination of oil extraction process and value of oil asset.

Experimental Methods for 2 & 3: Each sample was freshly prepared and tested (20μL sample volume) two times on the same day: 2 samples accumulating 4 scans per run using a CPMG pulse sequence with a scan duration of 4.0 seconds, a recycle delay of 2.5 seconds, and an echo period of 400 μs at ambient temperature. The data was fitted using a tri-exponential algorithm.

Objective 4: Analyze two sets of sulfur in diesel fuel ASI standards SFD7 (0 to 0.10%) and SFD10 (0 to 5%) to be able to grade diesel oil samples.

Objective 5: Analyze a petroleum distillate set of ONTA standards of Gasoline, Kerosene, and Diesel for differentiation.

Experimental Methods for 4 & 5: Each sample was freshly prepared and tested (50μL sample volume) two times on the same day using a CPMG pulse sequence with a scan duration of 2.0 seconds, a recycle delay of 2.5 seconds, and an echo period of 400 μs at ambient temperature. The data was fitted using mono and bi-exponential algorithms.

Conclusions

Using the WaveGuide Formula™ μNMR device WaveGuide was able to:

- Show that WaveGuide's μNMR can distinguish between synthetic motor oil grades and non-synthetic motor oils.
- Establish a linear slope for density and viscosity parameters in crude oils and differentiate motor oil samples from the original unused sample and from used samples at different time points.
- Classify wt.% asphaltene in crude oil samples in groups of high, medium and low concentration using T2a as a marker.
- Detect sulfur content in diesel fuel with a dynamic range of 50,000 ppm to 200 ppm.
- Measure distinct petroleum distillate fragments.



References

- Meiboom, S., and D. Gill. "Modified Spin-Echo Method for Measuring Nuclear Relaxation Times." *Review of Scientific Instruments*, vol. 29, no. 6, 1958, pp. 688-691. doi:10.1063/1.1713626.
- "Final Rule for Control of Air Pollution from Motor Vehicles; Tier 3 Motor Vehicle Emission and Fuel Standards." EPA-420-F-14-009. 2014. EPA, Environmental Protection Agency, 10 Sept. 2016. www.epa.gov/regulations-emissions-vehicles-and-engine/final-rule-control-air-pollution-motor-vehicles-tier-3.
- Belikov, Vladimir. "Shell Solar Split" uploaded by the Main Directorate of the Ministry of Internal Affairs of Russia for the Moscow Region, 18 December 2015. Retrieved from <https://www.youtube.com/watch?v=1aRG8h4750>
- Toma, Sebastian. "Engine Oil Counterfeiting Ring Uncovered in Russia." *Autoevolution*, 8 Jan. 2016, www.autoevolution.com/news/engine-oil-counterfeiting-ring-uncovered-in-russia-101064.html.
- Newsindianexpress.com/states/odisha/2016/oct/08/fake-engine-oil-manufacturing-unit-busted-in-odisha-1882591.html.
- Flammia, Dino. "Ni Puts the Brakes on Fraudulent Motor Oil." *New Jersey* 2015, New Jersey 1015.5 PM, 5 Dec. 2014, <http://nj1015.com/nj-puts-the-brakes-on-fraudulent-motor-oil/>.

Results

Table 1: T2 Values of 19 Commercially Available Motor Oils


SAMPLE	AVERAGE	STDEV	VISCOSITY
Bi-Exponential (T2A & T2B)	(4 samples)		mm ² /s @40°C
Pennzoil SAE 5W-30 Bottle #1 (Non-Synthetic)	T2A 0.03800 T2B 0.13665	0.00030 0.00018	63.4
Proline SAE 5W-30 Bottle #1 (Non-Synthetic)	T2A 0.04028 T2B 0.13555	0.00010 0.00026	67.1
Pennzoil SAE 10W-40 (Non-Synthetic)	T2A 0.02990 T2B 0.10598	0.00012 0.00030	104.7
Proline SAE 10W-40 (Non-Synthetic)	T2A 0.03038 T2B 0.10918	0.00004 0.00015	92.8
Castrol EDGE SAE 5W-30 GDI #1 (Synthetic)	T2A 0.04052 T2B 0.13633	0.00010 0.00022	62.64
Valvoline Advanced SAE 5W-30 #1 (Synthetic)	T2A 0.03632 T2B 0.12336	0.00009 0.00016	63.17
Mobil 1 5W-30 Bottle #1 (Synthetic)	T2A 0.03712 T2B 0.12530	0.00015 0.00036	63
Mobil 1 0W-30 Advanced Fuel Economy (Synthetic)	T2A 0.04708 T2B 0.15118	0.00011 0.00030	11.5
Castrol EDGE European SAE 0W-30 (Synthetic)	T2A 0.03598 T2B 0.13996	0.00008 0.00022	67.4
Castrol EDGE 10W-30 (Synthetic)	T2A 0.03292 T2B 0.11517	0.00006 0.00015	64.88
Mobil 1 5W-30 Bottle #2 (Synthetic)	T2A 0.03704 T2B 0.12525	0.00013 0.00037	63
Proline SAE 5W-30 Bottle #2 (Non-Synthetic)	T2A 0.04033 T2B 0.13600	0.00012 0.00036	67.1
Shell Rotella Gas Truck SAE 5W-30 (Synthetic)	T2A 0.04310 T2B 0.14061	0.00013 0.00041	66.4
Castrol EDGE SAE 5W-30 GDI #2 (Synthetic)	T2A 0.04062 T2B 0.137641	0.00010 0.00026	62.64
Castrol GTX Ultratec SAE 5W-30 (Synthetic)	T2A 0.03908 T2B 0.13093	0.00013 0.00032	61.34
Valvoline Advanced SAE 5W-30 #2 (Synthetic)	T2A 0.03822 T2B 0.12321	0.00010 0.00018	63.17
Castrol GTX Magnetic SAE 5W-30 (Synthetic)	T2A 0.03621 T2B 0.12545	0.00015 0.00045	61.77
Pennzoil Platinum SAE 5W-30 (Synthetic)	T2A 0.03803 T2B 0.12654	0.00006 0.00031	45.9
Pennzoil SAE 5W-30 Bottle #2 (Non-Synthetic)	T2A 0.03798 T2B 0.13082	0.00014 0.00028	63.4

Field Testing of Oil Samples




Table 2: 1/T2 vs Viscosity of Conditioned Motor Oils


Sample ID	Avg 1/T2 (s)	TIME (Hrs)	Viscosity (cSt)
WG-A1	35.40262 9.76348	0	90.96
WG-B1	35.0988 9.696218	0	90.86
WG-A2	32.39916 8.873823	1	81.71
WG-A3	30.08243 8.173641	16	67.83
WG-B2	32.3677 8.866191	1	82.03
WG-B3	29.8405 8.113327	8	67.94
WG-C1	32.23467 8.84228	1	82.06
WG-C2	29.9958 8.171938	16	68.05



Graph 1: T2 Comparison of Synthetic and Non-Synthetic Commercially Available Motor Oils
Comparing 1/T2 vs density, with 1/T2 on a log scale, one can see a correlation exists with the highest density materials displaying a fast T2, and with the less dense materials yielding a longer T2, as seen in Graph 3 (correlation coefficient of 0.94).




Graph 2: 1/T2 vs Viscosity of Conditioned Motor Oils



Graph 4: T2 of Crude Oil vs. wt% Asphaltene



SAMPLE	AVERAGE	STDEV	density @ 15°C(g/ml)	Asphaltene (wt%)	Viscosity (kinematic) 40°C	Viscosity (absolute) 80°C
Tri-Exponential (T2A, T2B & T2C)	(2 samples)					
WG1	T2A 0.00506 T2B 0.02178 T2C 0.06882	0.00006 0.00003 0.00056	0.9478	0.05	101.0	15.2
WG2	T2A 0.00567 T2B 0.12619 T2C 0.41114	0.00016 0.00097 0.00614	0.8779	0.07	10.8	2.5
WG3	T2A 0.00268 T2B 0.01132 T2C 0.03158	0.00004 0.00012 0.00011	0.9125	0.55	37.9	8.6
WG4	T2A 0.00246 T2B 0.12296 T2C 0.38697	0.00034 0.00134 0.00333	0.8786	0.45	11.9	3.7
WG5	T2A 0.00283 T2B 0.01005 T2C 0.02669	0.00005 0.00022 0.00033	0.9422	0.92	219.2	27.0
WG6	T2A 0.01843 T2B 0.09642 T2C 0.33060	0.00001 0.00029 0.00213	0.8871	0.05	21.4	5.2
WG7	T2A 0.00882 T2B 0.04074 T2C 0.12884	0.00002 0.00014 0.00006	0.9123	0.85	35.5	8.4
WG8	T2A 0.00906 T2B 0.04233 T2C 0.11953	0.00008 0.00014 0.00015	0.8826	1.60	25.3	5.6
WG9	T2A 0.01593 T2B 0.08294 T2C 0.26809	0.00011 0.00059 0.00175	0.8846	0.84	17.6	3.8
WG10	T2A 0.00517 T2B 0.00141 T2C 0.06283	0.00012 0.00051 0.00134	0.9341	2.10	79.3	14.6
WG11	T2A 0.01525 T2B 0.07345 T2C 0.21109	0.00029 0.00139 0.00902	0.8751	1.60	9.9	3.1
WG12	T2A 0.01440 T2B 0.06946 T2C 0.19461	0.00005 0.00064 0.00147	0.8758	1.76	10.3	3.2
WG13	T2A 0.09123 T2B 0.43143 T2C 0.25217	0.02497 0.10619 0.15257	0.8218	0.32	2.7	1.1
WG14	T2A 0.14738 T2B 0.52805 T2C 0.12620	0.00438 0.1079 0.2019	0.8429	0.00	2.6	1.1
WG15	T2A 0.15461 T2B 0.73061 T2C 0.17352	0.00321 0.00207 0.00298	0.8531	0.19	1.6	0.8
WG16	T2A 0.08611 T2B 0.40263 T2C 0.24248	0.00252 0.09429 0.16123	0.8204	0.03	3.2	1.2
WG17	T2A 0.06638 T2B 0.19047 T2C 0.14777	0.00015 0.00017 0.00022	0.9041	0.27	23.5	5.6
WG18	T2A 0.05267 T2B 0.13459 T2C 0.02168	0.00334 0.1009 0.00037	0.8086	0.02	2.3	1.0
WG19	T2A 0.10626 T2B 0.15461 T2C 0.38564	0.00184 0.00321 0.00453	0.8917	0.69	14.5	4.2
WG20	T2A 0.00305 T2B 0.01099 T2C 0.02987	0.00014 0.00038 0.00037	0.9366	6.68	79.6	15.6
WG21	T2A 0.00501 T2B 0.13147 T2C 0.36188	0.00036 0.00264 0.00882	0.8430	2.18	4.7	1.9
WG22	T2A 0.05508 T2B 0.26660 T2C 0.17746	0.00087 0.00334 0.00651	0.8369	0.34	4.0	1.5
WG23	T2A 0.03660 T2B 0.13700 T2C 0.34780	0.00081 0.00511 0.00835	3.33	5.2	1.7	
WG24	T2A 0.02373 T2B 0.11127 T2C 0.28565	0.00045 0.00163 0.00265	0.8531	1.81	5.6	2.1
WG25	T2A 0.00344 T2B 0.01375 T2C 0.03879	0.00069 0.00429 0.00298	0.9309	7.87	88.1	19.2
WG26	T2A 0.00448 T2B 0.11213 T2C 0.00394	0.00050 0.00085 0.00009	0.8730	2.76	14.5	4.0
WG27	T2A 0.01628 T2B 0.05418 T2C 0.03141	0.00021 0.00029 0.00031	0.9311	4.96	67.9	14.5
WG28	T2A 0.16715 T2B 0.47840 T2C 0.00353	0.00247 0.00721 0.00004	0.8572	1.65	6.4	2.3
WG29	T2A 0.01422 T2B 0.00380 T2C 0.00380	0.00018 0.00003 0.00003	0.9354	8.13	125.2	23.1
WG30	T2A 0.00647 T2B 0.03598 T2C 0.02940	0.00724 0.00297 0.00011	0.9334	7.05	71.0	17.7
WG31	T2A 0.15356 T2B 0.45768 T2C 0.00848	0.00203 0.00848 0.00048	0.8597	1.74	5.7	2.4

Graph 6: T2 of Crude Oil vs. Absolute Viscosity




Graph 5: T2 of Crude Oil vs. Kinematic Viscosity

Results Viscosity: Comparing 1/T2 vs kinematic viscosity (the time it takes for a specific volume of oil to flow through a capillary tube) one can see a linear-like slope exists (correlation coefficient of 0.94) with the more viscous oil samples displaying a fast T2A, and with less viscous materials yielding a longer T, (seen in Graph 5). Comparing 1/T2 vs absolute viscosity (the force needed by a fluid to overcome its own internal molecular friction so that it can flow) one can see a linear-like slope exists (correlation coefficient of 0.95) with the more viscous oil samples displaying a fast T2, and the less viscous materials yielding a longer T2 (seen in Graph 6).




Graph 7: T2 of ASI Diesel Oil vs. Concentration

Results ASI Diesel Oil: Data was analyzed and plotted using the average T2 measurements of samples collected in triplicate with a single, and a bi-exponential fit for the ASI standards. The WaveGuide Formula™ was able to effectively detect in a linear fashion the known set of sulfur standards with only a single exponential fit shown (Graphs 7-8). The dynamic limit of detection using the small molecule sulfur gas chromatography ASI standard was 500,000 ppm to 200 ppm of total sulfur.



Graph 8: T2 of ASI Diesel Oil vs. Concentration



Graph 9: T2 of Petroleum Distillates

Results ONTA Petroleum Distillate Standards: The data was analyzed using average T2 measurements of samples tested in triplicate with single, and bi-exponential fits for the ONTA standards. The WaveGuide Formula™ was able to effectively detect the difference in the distillates standards with a single exponential shown (Graph 9). Gasoline has the longest T2 and the crude diesel sample had the shortest T2.